EN DeVilbiss® SmartLink™ Module Instruction Guide, Model DV5M

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TABLE OF CONTENTS

Introduction.....	3
Symbol Definitions	3
Important Safeguards	4
Set Up	5
Important Parts of Your SmartLink Module.....	5
Connection to the CPAP	5
Optional Oximeter Connection	6
Operating Instructions.....	8
Navigating Menus and Onboard Reporting	8
Data Retrieval	10
Accessories.....	12
Cleaning Instructions	12
Service Instructions	12
Troubleshooting.....	13
Specifications	14
Declaration of Conformity.....	15

INTRODUCTION

The DeVilbiss SmartLink Module is an optional accessory to the DV5X series CPAP devices. This SmartLink device was supplied to you by your homecare provider as part of your overall positive airway pressure therapy for obstructive sleep apnea (OSA). SmartLink is a helpful tool for your physician and/or homecare provider because it allows them to see the effectiveness of the therapy and how frequently and consistently you use the device.

The information SmartLink collects helps bring to light potential problems or challenges you are facing that your clinician might be able to help you with in getting used to therapy.

SmartLink is an advanced technology module that works together with your CPAP (Continuous Positive Airway Pressure) device.

This system is unique in healthcare because many times a clinician will prescribe a therapy and send you home with it and never really know how effective it is. SmartLink is a great benefit to you because it helps your clinician keep track of what's going on with your condition without you having to make extra trips to the clinician's office.

INDICATIONS FOR USE

The DeVilbiss SmartLink System can only be used in conjunction with the DV51 and DV54 Series CPAP Systems for follow up of obstructive sleep apnea in patients weighing above 30 kg on nasal CPAP therapy.

SYMBOL DEFINITIONS

Attention - Consult Accompanying Documents

Type BF equipment— applied part

Class II electrical protection— double insulated

Data Port Input/Output

Memory Card Insertion Orientation

Oximeter Connection Port

This device contains electrical and/or electronic equipment that must be recycled per EC Directive 2002/96/EC- Waste Electrical and Electronic Equipment (WEEE)

IMPORTANT SAFEGUARDS

Please read all instructions carefully before use. Important information is highlighted by the following terms:

DANGER	Urgent safety information for hazards that will cause serious injury or death.
WARNING	Important safety information for hazards that might cause serious injury.
CAUTION	Information for preventing damage to the product.
NOTE	Information to which you should pay special attention.

DANGER

- **Electric Shock Hazard-Do not use while bathing.**
- **Electric Shock Hazard-Do not immerse this device into water or any other liquid.**
- **Electric Shock Hazard-Do not attempt to open or remove the enclosure.** There are no user-serviceable components inside. Contact your home care provider if your device needs servicing. Opening or tampering with the product will void the warranty.

WARNING

- **Use the SmartLink Module only for its intended use as described in this manual.**
- **Be sure to read and understand all safety instructions supplied with your SmartLink Module device.**
- **Use only accessories recommended by DeVilbiss.**
- **Use only data cards provided by your homecare provider. Other Secure Digital (SD) cards may not be compatible with the system or may not perform properly.**
- **Do not expose the device to water.**
- **If water has been spilled onto the SmartLink Module or it has been submerged into water, unplug the power cord from the power source. Allow the device to completely dry before use.**
- **If the device has been dropped, refer to the troubleshooting guide for instructions. Contact your home care provider for evaluation of the device.**
- **To avoid electric shock, always unplug power cord from wall outlet power**

source when performing cleaning.

- Do not insert any foreign objects into any openings.

CAUTION—The circular data port and oximeter connectors located on the side of the SmartLink Module are used to attach accessories to the device. The connectors must only be used with accessories approved for use by DeVilbiss. Do not attempt to attach any other device to these connectors as it may damage the CPAP, SmartLink Module, and/or the device you are attempting to attach.


Be sure that you read and understand the information in this instruction guide. The information in this guide will help you understand the operation of the Smartlink Module and how to properly maintain it. By following these instructions and the advice of your physician, your SmartLink Module will become an effective addition to your CPAP therapy.

SAVE THESE INSTRUCTIONS

SET UP


Important Parts Of Your SmartLink Module

1. Data card slot
2. Release button
3. Secondary accessory data connection port
4. Oximeter connection port
5. Primary data connector
6. Locking tabs (3 tabs)
7. Instruction Guide (not shown)


Connection To The CPAP

NOTE—Your CPAP may already have been provided to you with the SmartLink attached. If you have received your module separately from the CPAP, follow the steps below to properly attach the module to the CPAP.

1. Locate the primary data connector on the back of the SmartLink Module. Align and engage the primary data connector into the circular DIN connector on the back of the CPAP. At the same time, align and engage the three locking tabs on the module with the three locking recesses on the back surface of the CPAP **Figure 1**.
2. **Firmly** press the module onto the back of the CPAP until the module is fully secured and locked into position **Figure 2**. An audible click may be heard confirming that the module is properly attached to the CPAP. Visually confirm that there is not a gap between the module and the CPAP.

Removal From The CPAP

To remove the module from the CPAP:

1. Fully depress the release button on the side of the module.
2. Pull the module straight back and away from the CPAP.

CAUTION-Do not pull up on the module as this may damage the data connector.

Optional Oximeter Connection

The SmartLink Module provides oximeter connection capability to allow an approved oximeter to be used with the CPAP/module system. Connect the oximeter cable into the oximeter connection port located on the side of the module. If you are using an oximeter, please refer to the assembly and operation instructions that came with the oximeter **Figure 3**.

You can see your oximeter readings by using the SmartLink menu (see SmartLink Menu Structure) at any time.

Message shown when the oximeter is first connected and logging starts. This message will remain on display for five seconds. The pulse oximeter menu (SpO_2 /PR) will then be displayed.

Oximeter
Connected

Message shown when the oximeter is removed. This message will remain on display for five seconds. Multiple notifications will occur until logging stops.

Oximeter
Removed

SmartLink Menu Structure


To display or remove SmartLink menu, press Left and Right Arrow Keys together. Use Arrows to navigate through the menu structure as shown (will follow dotted line).

Navigation:

- ▶ Right Arrow - navigates right through menu items.
- ◀ Left Arrow - navigates left through menu items.
- ▲ Up Arrow - navigates up through menu items.
- ▼ Down Arrow - navigates down through menu items.

To Exit SmartLink menu (choose one):

1. Turn blower on or off
2. Press Up Arrow Key at top level Quick View Menu
3. Press Left and Right Arrow Keys together.

OPERATING INSTRUCTIONS

Navigating Menus and Onboard Reporting

While the SmartLink Module is attached, it will collect important data about how the device is being used and what kind of therapy it is delivering for you. Much of the information that is collected is available to you through the onboard menu on your CPAP display. The previous table shows how to navigate through the menus.

When the SmartLink is connected correctly, the display on the CPAP will show this message. **NOTE**—if a data card is inserted in the module before attaching the module to the CPAP, then card detection messages may be displayed.

SmartLink Module
Ready for Use

To enter the SmartLink menu, press the LEFT and RIGHT Arrow keys on the CPAP keypad at the same time. To exit the SmartLink menu, press the same key sequence again or simply turn the CPAP off and on.

SmartLink Menu
Press ◀▶ To Exit

The next menu item is the “Quick View” submenu. To enter the Quick View submenus, press the DOWN Arrow key. **NOTE**—pressing the UP Arrow key at the Quick View menu will exit the SmartLink menu.

Quick View Menu
Press ▼

Quick View Submenu Navigation:

Navigate through the Quick View submenu by using the UP or DOWN Arrow keys to change between 1, 7, or 30 day averages and RIGHT or LEFT Arrow keys to cycle through the items within the selected time period.

The average pressure provided during the selected time period (1, 7, or 30 days).

Average Pressure
30 day avg 10.1

The Apnea/Hypopnea Index (AHI) is the number of apneas and hypopneas detected per hour. The average AHI is displayed for the time period selected (1, 7, or 30 days).

AHI
30 day avg: 2.0

The Leak Percentage is the percentage of therapy time that the system leak was above 95 Liters/Minute (L/Min). A system leak above 95 L/Min, in most cases, is considered excessive and is usually an indication of a poor mask fit that should be corrected

Leak %
30 day avg: 5.0

to ensure optimal therapy. The average leak percentage is displayed for the time period selected (1, 7, or 30 days).

The Non-Responding Event Index (NRI) is the number of non-responding events detected per hour. NRI is a clinical measurement that may be monitored by your clinician. The average NRI is displayed for the time period selected (1, 7, or 30 days).

NRI
30 day avg: 1.0

The Exhale Puff Index (EPI) is the number of exhale puff events detected per hour. Exhale puffing is caused by inhaling through the patient circuit and exhaling through the mouth. EPI is a clinical measurement that may be monitored by your clinician. The average AHI is displayed for the time period selected (1, 7, or 30 days).

EPI
30 day avg: 1.0

The average number of hours of therapy provided during the selected time period (1, 7, or 30 days).

Hours Per Day
30 day avg: 8.0

The Adherence Percentage is the percentage of days of usage in which the system was used for a total of 4 or more hours per day. The adherence percentage is displayed for the time period selected (1, 7, or 30 days).

Adherence %
100% (1234/1234)

Press the UP Arrow key on the CPAP keypad to exit the Quick View submenu and return to the SmartLink menu structure.

Exit Menu
Press ▲

Press the UP or DOWN arrow key on the CPAP keypad to adjust the volume for the audible alert notification.

Volume
100% (Adjust ▲▼)

This menu displays the serial number of the SmartLink module connected to the CPAP.

Serial Number
CM123456

This menu displays the operating firmware version of the SmartLink module connected to the CPAP.

Firmware Version
V0.06 25/04/2008

This menu displays the pulse oximetry (SpO₂) and pulse rate (PR) information when a pulse oximeter is used with the system. If a pulse oximeter is not connected to the module “***” will be displayed for the SpO₂ and PR values.

SpO₂: 100
PR: 101

Data Retrieval

Your clinician may have a standard protocol to follow-up on your therapy to find out how the device is working and to identify any problems or issues you may be having. Your home care provider will mail a SmartLink Data Card to you at certain intervals and ask you to download the data onto the card and mail it back. It is a very simple process:

1. Ensure that the SmartLink module is properly connected to the CPAP and that the CPAP is connected to an appropriate power source.
2. Orient that card as shown in **Figure 4**. Insert and press the card into the card slot located on the top of the module until the card clicks into place **Figure 5**.
3. Ensure that the card is detected (view message or audible alert) and allow the data to transfer to the card. Depending on the amount of data, this could take several minutes. Do not remove the card from the module until complete. Observe the message(s) on the display to confirm that the transfer was successful (see Messages below). If necessary, contact your provider and report any errors that may have occurred.
4. If instructed by your provider, eject the card from the card slot by pressing down on the card until it is released and remove the card from the slot.
5. Place the card in the provided mailer and return it to the provider.

Messages

The following messages could be displayed on the CPAP during the data retrieval process:

Message shown when a card is first inserted into the SmartLink Module and data is being transferred. This message will remain on the display until transfer is completed.

Card Detected
Please Wait . . .

Message shown when a card is first inserted into SmartLink Module and there is no new data to transfer. Will remain on display for ten seconds.

Card Detected
Ready For Use

Message shown after completion of card data transfer. This message will remain on the display for ten seconds.

Card Transfer
Successful

Message shown if an error is encountered during card data transfer. This message will remain on the display until acknowledged by pressing any key on the keypad or a 15 minute timeout occurs.

Card Error
Contact Provider

Message shown if an error is encountered while updating device settings. This message will remain on the display until acknowledged by pressing any key on the keypad or a 15 minute timeout occurs. **NOTE**—*If this message is displayed, contact your clinician.*

Setting Error
Contact Provider

Message shown after removing the card from the module. This message will remain on the display for ten seconds.

Card Removed
Send to Provider

Message shown when a new module with mismatched data is detected. This message will remain on the display until the module is removed. Contact your provider to notify them of the situation and for instructions on how to resolve the issue. If this situation is not resolved, the module will not be able to record usage and therapy data. However, although this message will remain on the display with the module attached, the CPAP will still be functional and provide therapy.

Setup Required
Contact Provider

Other Messages

In the unlikely event that a fault should occur, this message will be shown when a module fault condition is encountered (Clock Error, Memory Error, Setting Error, etc).

E00 Module Fault
Contact Provider

Message shown when firmware is being updated. This message will remain on the display until the update is complete.

Firmware Update
Please Wait . . .

Occasionally a firmware update may be provided for your device via SmartLink Data Card. The firmware update process is automatic and may take a few seconds to complete. Once the update is complete, the system will return to normal operation.

ACCESSORIES

Oximeter Kit	DV5M-697
Oximeter Finger Clip Sensor	8000AA
Oximeter Ear Clip Sensor	8000Q
nVision® Software.....	9054D-699

CLEANING INSTRUCTIONS

WARNING

To avoid electric shock, always unplug power cord from the wall outlet power source.

CAUTION—*Never rinse or place the device in water. Never allow liquids to get into or around any of the ports or openings; doing so will result in product damage. If this occurs, do not use the device. Contact your provider for service.*

NOTE—*The module should be cleaned as needed.*

1. Turn off the CPAP and unplug it from the power source.
2. Use a soft damp cloth to clean the surfaces of the device.
3. Allow the device to dry completely before plugging it into a power source.

SERVICE INSTRUCTIONS

The SmartLink Module does not require regular servicing.

TROUBLESHOOTING

Problem	Probable Cause	Remedy
Module is not recognized when connected to the CPAP.	1. The CPAP is not connected to a power source.	1. Connect the CPAP to an appropriate power source and verify that the module is recognized.
	2. Module is not fully attached to the CPAP.	2. Remove the module from the CPAP, realign the primary data connector with the circular DIN connector of CPAP, and press firmly until the module is locked into position. Be sure there is not a gap between the module and the CPAP.
Card not recognized when inserted into the module card slot.	1. Card not oriented properly.	1. Refer to the symbol on the top of the module or Figure 4 in this guide for proper card orientation.
	2. Card not fully inserted in the module card slot.	2. Press the card fully into the card slot until the card clicks into position.
	3. Module is not fully attached to the CPAP.	3. Remove the module from the CPAP, realign the primary data connector with the circular DIN connector of CPAP, and press firmly until the module is locked into position. Be sure there is not a gap between the module and the CPAP.
Card error.	1. Wrong card type.	1. Only use a card supplied by your provider. The card must be an SD card with FAT16 format.
	2. Card removed during data transfer.	2. Card must not be removed during data transfer. Data transfer interruption may corrupt the card and prevent the card from working correctly. Reinsert card and allow data transfer to continue until completed.
	3. Card corrupt or faulty.	3. Return card to provider for replacement.

Problem	Probable Cause	Remedy
Settings error.	1. Card removed during data transfer.	1. Card must not be removed during data transfer. Data transfer interruption may corrupt the card and prevent the card from working correctly. Reinsert card and allow data transfer to continue until completed. If this error is not corrected, contact provider for resolution.
	2. Card corrupt or faulty.	2. Return card to provider for replacement. If this error is not corrected, contact provider for resolution.
Setup required.	1. Module not set up properly.	1. Notify provider of the situation and follow provider instructions.
Oximeter data (SpO ₂ and PR) not displayed.	1. Oximeter not properly connected to the SmartLink module.	1. Verify connection between the oximeter connector and module oximeter connection port.
	2. Oximeter probe not properly attached to user.	2. Reference the operating instructions that came with oximeter.
CPAP is not fully functional.	1. Module may be defective.	1. Remove the module from the CPAP and verify that CPAP is fully functional. If it functions properly without the module connected, the module may be defective. Contact provider. If it still is not fully functional, consult the CPAP instruction guide for troubleshooting.

DEVILBISS SMARTLINK SPECIFICATIONS

Size: 3.9" H x 3.1" W x 1.4" D (9.9cm H x 7.9cm W x 3.6cm D)

Weight: 0.30 lbs. (0.14 Kg)

Electrical Rating:

Maximum Power Consumption (SmartLink Module only): 4 W

Voltage and Current: Power supplied by the DV5X Series Device

Class II Equipment; Type BF Applied Parts; Continuous Operation

Operating Conditions:

Temperature Range: 41°F to 104°F (5°C to 40°C)

Humidity Range: 0 to 95% R.H. non-condensing

Atmospheric Conditions: Sea level to 9,000' (2,743 M)

Transportation and Storage Conditions:

Temperature Range: -40°F to 158°F (-40°C to +70°C)

Humidity Range: 0 to 95% R.H. non-condensing

DECLARATION OF CONFORMITY

Manufacturer: Sunrise Medical dba DeVilbiss Healthcare

Address: 100 DeVilbiss Drive, Somerset, PA 15501

Product Designation: CPAP Accessory

Type, Model: DV5M SmartLink Module

For DeVilbiss IntelliPAP or SleepCube CPAP

We herewith declare that the above mentioned product complies with the requirements of EC Directive 93/42/EEC and the following:

Class: IIa, Rule 2

Quality System
Standards Applied: ISO 13485:2003

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DeVilbiss Healthcare Pty. Limited
DeVilbiss Healthcare SAS
DeVilbiss Healthcare GmbH
DeVilbiss Healthcare LLC • 100 DeVilbiss Drive • Somerset, PA 15501 • USA